


BRUSHES

GASTROINTESTINAL PRODUCTS

Microbiology Brushes: Protected Specimen Brush

- Teflon sheathing for ease of insertion and manipulation during procedures
- Spring-loaded handle and permanent tear drop shaped tip ensures a protective seal
- Latch hook for locking handle in open brush position
- Ergonomic handle provides operator comfort
- Brush configuration ensures a high specimen yield
- Packaged sterile for single use: 5 per box



MICROBIOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4352	2.0mm	1.7mm	220cm	Green Swirl

PULMONARY PRODUCTS

MICROBIOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4310	1.0mm	1.7mm	100cm	Blue Swirl
4320	2.0mm	1.7mm	100cm	Blue Swirl





Modified Saccomanno™ Cell Transport System

- Allows for quick washing of brushes, forceps and needles
- Immediate cell fixation
- Convenient self-contained vial with patient ID tags
- Packaged 10 per box

MODIFIED SACCOMANNO™ CELL TRANSPORT SYSTEM

Catalog Number	Description
4101	Cell Transport System

PULMONARY PRODUCTS



Nasal Cytology Brush

- Tough stainless steel, high tempered wire to prevent kinking
- Bullet tip prevents mucosal trauma
- Color-coded handle provides quick and easy visual identification
- Packaged sterile for single use: 10 per box

NASAL CYTOLOGY BRUSH

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4290	2.0mm	N/A	10cm	Blue Swirl
4291	3.0mm	N/A	10cm	Blue Swirl



GASTROINTESTINAL PRODUCTS

Cytology Brushes

- Teflon sheathing for ease of insertion and manipulation during procedures
- High tempered stainless steel wire to prevent kinking
- Radiopaque ball tip prevents mucosal trauma
- Thumb ring for easy operation
- Endoscope compatible lengths
- Packaged sterile for single use: 10 per box



CYTOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4102	3.0mm	1.7mm	160cm	Green
4104	3.0mm	1.7mm	200cm	Green
4106	3.0mm	1.7mm	230cm	Green
4108	3.0mm	2.3mm	160cm	Green
4110	3.0mm	2.3mm	200cm	Green
4112	3.0mm	2.3mm	260cm	Green

PULMONARY PRODUCTS

CYTOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4202	1.0mm	1.7mm	120cm	Blue
4204	2.0mm	1.7mm	120cm	Blue
4206	3.0mm	1.7mm	120cm	Blue





Cleaning Brushes: Disposable

GOOD-BI-BRUSH (2-BRUSH CLEANING SYSTEM)

- Unique tapered 2-brush system provides dual cleaning action with one pass
- Blunt, looped tip brush protects delicate channels
- Improved sheath flexibility for easy passage through working channel and umbilical exit port
- Vivid color cannula offers easy identification in cleaning solution
- Endoscope compatible lengths
- Packaged for single use: 25 per box

GOOD-BI-BRUSH (2-BRUSH CLEANING SYSTEM): DISPOSABLE

Catalog Number	Brush Taper	Brush Length	Sheath	Length	Working Channel
4660	3-5mm x 4-6mm	2cm x 1cm	1.7mm	260cm	2.0mm & Larger

Single Cleaning Brush System

- Unique tapered brush provides greater channel cleaning

SINGLE CLEANING BRUSH SYSTEM: DISPOSABLE

Catalog Number	Brush Taper	Brush Length	Sheath	Length	Working Channel
4676	3-5mm	2cm	1.7mm	260cm	2.0mm & Larger
4677	3-5mm	1cm	1.7mm	260cm	2.0mm & Larger

PULMONARY PRODUCTS

- Unique tapered 2-brush system provides dual cleaning action with one pass

GOOD-BI-BRUSH (2-BRUSH CLEANING SYSTEM): DISPOSABLE

Catalog Number	Brush Taper	Brush Length	Sheath	Length	Working Channel
4665	3-5mm x 4-6mm	2cm x 1cm	1.7mm	100cm	2.0mm & Larger





HOBBS MEDICAL, INC.

CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

Orders & Returns

Product Information
Pricing
Specification Changes