

GASTROINTESTINAL PRODUCTS

Microbiology Brushes: Protected Specimen Brush

- Teflon sheathing for ease of insertion and manipulation during procedures
- Spring-loaded handle and permanent tear drop shaped tip ensures a protective seal
- Latch hook for locking handle in open brush position
- Ergonomic handle provides operator comfort
- Brush configuration ensures a high specimen yield
- Packaged sterile for single use: 5 per box



MICROBIOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4352	2.0mm	1.7mm	220cm	Green Swirl

PULMONARY PRODUCTS

MICROBIOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4310	1.0mm	1.7mm	100cm	Blue Swirl
4320	2.0mm	1.7mm	100cm	Blue Swirl

