

Biopsy Forceps – Single Use



Oval Fenestrated



Oval Fenestrated with Spike

Gastroscopy & Colonoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0222-B	Oval Fenestrated	2.3	160	2.8
NE0222-C	Oval Fenestrated	2.3	180	2.8
NE0222D-B	Oval Fenestrated with Spike	2.3	160	2.8
NE0222D-C	Oval Fenestrated with Spike	2.3	180	2.8
NE0222D-G	Oval Fenestrated with Spike	2.3	230	2.8
NE0222-G	Oval Fenestrated	2.3	230	2.8
NE0230D-G	Oval Fenestrated with Spike	3.0	230	3.3
NE0230-G	Oval Fenestrated	3.0	230	3.3
NE1618-B	Oval Fenestrated	1.8	160	2.0
NE1618-C	Oval Fenestrated	1.8	180	2.0



Alligator



Alligator with Spike

Gastroscopy & Colonoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0422D-G	Alligator with Spike	2.3	230	2.8
NE0422-G	Alligator	2.3	230	2.8
NE0430-G	Alligator	3.0	230	3.3



LEF-GBF-SU-00-EN / LL-0359 Vers. 2.0 (10/2020)

Biopsy Forceps

Single Use



Standard and Serrated (Alligator) Cups

Precise and clean tissue acquisition

Versions with Teflon Coated Metal Spiral

Friction minimized passage through accessory channel



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
 Email: sales@endo-flex.de • www.endo-flex.de

Biopsy Forceps – Single Use

Cup Shapes



Bronchoscopy

Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0418-A	Alligator	1.8	120	2.0
NE1618-A	Oval Fenestrated	1.8	120	2.0



Enteroscopy


Uncoated stainless steel spiral

Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE0222-M	Oval Fenestrated	2.3	280	2.8
NE1618-M	Oval Fenestrated	1.8	280	2.0



Biopsy Forceps – Single Use

Teflon Coated Stainless Steel Spiral



Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NEX0222-B	Oval Fenestrated	2.3	160	2.8
NEX0222D-B	Oval Fenestrated with Spike	2.3	160	2.8
NEX0222D-G	Oval Fenestrated with Spike	2.3	230	2.8
NEX0222-G	Oval Fenestrated	2.3	230	2.8
NEX0230-G	Oval Fenestrated	3.0	230	3.3



Teflon Coated Stainless Steel Spiral



Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NEX0422-G	Alligator	2.3	230	2.8
NEX0430-G	Alligator	3.0	230	3.3



Hot Biopsy Forceps – Single Use



Hot Biopsy Forceps – Teflon Coated Stainless Steel Spiral



Catalogue No.	Cup Shape	Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE6122-G	Oval Fenestrated	2.3	230	2.8

