

Cytology Brushes

Single Use



Brush with Metal Tip

Efficient tube closure to secure smear yield

Double Lumen Version

Compatible with guide wire

Cytology Brushes – Single Use



Cytology Brush – Single Lumen

Catalogue No.	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
E4218-A	120	1.8	2.0
E4218-C	180	1.8	2.0
E4219-G	230	1.9	2.2
E4222-G	230	2.3	2.8



Cytology Brush – Single Lumen

Catalogue No.	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
NE4218-A	120	1.8	2.0
NE4218-C	180	1.8	2.0
NE4222-G	230	2.3	2.8



Cytology Brush – Double Lumen

Compatible with guide wire .035"

Catalogue No.	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
E4222DL-C	180	2.7	3.2



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
 Email: sales@endo-flex.de • www.endo-flex.de