

Injection Needles

Single Use



LEF-GSN-SU-01-EN / LL-0356 Vers. 1.0 (08/2020)

Metal Tip Protection

Secured position of needle tip in the tube and safe needle guidance in angulated endoscope positions

Transparent Needle Tube

Flow control of therapeutics and fluids

Needle Lock Mechanism


Controlled puncture in targeted tissue and vessels

ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de



Injection Needles – Single Use



Injection Needle STANDARD-Series

Removable safety clip; Luer-Lock-mechanism

Catalogue No.	Needle Diameter (mm / Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
ET2518-C4	0.6/23	4	1.8	180	2.0
ET2518-C6	0.6/23	6	1.8	180	2.0
ET2518-G4	0.6/23	4	1.8	230	2.0
ET2522-B4	0.7/22	4	2.3	160	2.8
ET2522-B6	0.7/22	6	2.3	160	2.8
ET2522-C4	0.7/22	4	2.3	180	2.8
ET2522-C6	0.7/22	6	2.3	180	2.8
ET2522-G4	0.7/22	4	2.3	230	2.8
ET2522-G6	0.7/22	6	2.3	230	2.8
ET2522-M6	0.7/22	6	2.3	280	2.8
ET2522-O6	0.7/22	6	2.3	300	2.8
ET2622-B6	0.9/20	6	2.3	160	2.8
ET2622-C6	0.9/20	6	2.3	180	2.8
ET2622-G6	0.9/20	6	2.3	230	2.8



Injection Needle PRO-Series

Metal tip protection; needle lock mechanism (one hand)

Catalogue No.	Needle Diameter (mm / Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NET2422-C4	0.5/25	4	2.3	180	2.8
NET2422-G4	0.5/25	4	2.3	230	2.8
NET2422-G6	0.5/25	6	2.3	230	2.8
NET2518-C4	0.5/25	4	1.8	180	2.0
NET2522-B4	0.7/22	4	2.3	160	2.8
NET2522-B6	0.7/22	6	2.3	160	2.8
NET2522-C4	0.7/22	4	2.3	180	2.8
NET2522-C5	0.7/22	5	2.3	180	2.8
NET2522-C6	0.7/22	6	2.3	180	2.8
NET2522-G4	0.7/22	4	2.3	230	2.8
NET2522-G5	0.7/22	5	2.3	230	2.8
NET2522-G6	0.7/22	6	2.3	230	2.8
NET2622-B6	1.0/19	6	2.3	160	2.8
NET2622-C4	1.0/19	4	2.3	180	2.8
NET2622-C6	1.0/19	6	2.3	180	2.8
NET2622-G4	1.0/19	4	2.3	230	2.8
NET2622-G5	1.0/19	5	2.3	230	2.8
NET2622-G6	1.0/19	6	2.3	230	2.8

