


Polypectomy / Mucosectomy Snares

Single Use



LEF-GPS-SU-00-EN / LL-0363 Vers. 1.0 (08/2020)



Hexagonal Loop



Polypectomy Snare

Hexagonal loop, braided wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE372214-G	Hexagonal	15	2.35	2.8	230
NOE372216-G	Hexagonal	25	2.35	2.8	230
NOE372217-G	Hexagonal	35	2.35	2.8	230



Asymmetrical Loop



Polypectomy Snare

Asymmetrical loop, braided wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE3222-G	Asymmetrical	15	2.35	2.8	230
NOE3322-G	Asymmetrical	25	2.35	2.8	230



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de

Versions with 360° Rotation
Easier polyp snaring

Monofilament Loop Versions
Precise and smooth cutting

Rotatable Polypectomy Snares – Single Use




Oval Loop

Polypectomy Snare

Oval loop, braided wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE341816-C	Oval	25	1.8	2.0	180
NOE342213-G	Oval	10	2.35	2.8	230
NOE342214-C	Oval	15	2.35	2.8	180
NOE342214-G	Oval	15	2.35	2.8	230
NOE342214-M	Oval	15	2.35	2.8	280
NOE342216-C	Oval	25	2.35	2.8	180
NOE342216-G	Oval	25	2.35	2.8	230
NOE342217-C	Oval	35	2.35	2.8	180
NOE342217-G	Oval	35	2.35	2.8	230
NOE342217-M	Oval	35	2.35	2.8	280



Oval Loop

Polypectomy Snare

Oval loop, monofilament wire, rotatable

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
NOE342212M-G	Oval	6	2.35	2.8	230
NOE342213M-G	Oval	10	2.35	2.8	230
NOE342214M-C	Oval	15	2.35	2.8	180
NOE342214M-G	Oval	15	2.35	2.8	230
NOE342216M-C	Oval	25	2.35	2.8	180
NOE342216M-G	Oval	25	2.35	2.8	230
NOE342217M-G	Oval	35	2.35	2.8	230



Polypectomy Snares – Single Use



Oval Loop

Polypectomy Snare

Oval loop, braided wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE342215-G	Oval	20	2.3	2.8	230
OE342216-G	Oval	25	2.3	2.8	230
OE342217-C	Oval	35	2.3	2.8	180
OE342217-G	Oval	35	2.3	2.8	230



Oval Loop

Polypectomy Snare

Oval loop, monofilament wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE342214M-G	Oval	15	2.3	2.8	230



Asymmetrical Loop

Polypectomy Snare

Asymmetrical loop, braided wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE3322-G	Asymmetrical	25	2.3	2.8	230



Mucosectomy Snares – Single Use



Oval Loop

Mucosectomy Snare

Oval loop, braided wire

Catalogue No.	Loop Shape	Loop Diameter (mm)	Sheath Diameter (mm)	Min. Accessory Channel (mm)	Working Length (cm)
OE342246-G	Oval	25	2.3	2.8	230

