

Retrieval Devices

Single Use



LEF-GFB-SU-01-EN / LL-0362 Vers. 1.0 (07/2020)

Polyp and Foreign Body Retrievers

Strong & wide opening graspers
Large variety of applications

Foreign Body Forceps

Coated metal spiral
Friction minimized passage through accessory channel



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de


Retrieval Devices – Single Use

Foreign Body Forceps – Single Use



Foreign Body Retriever

Teflon tube




3 Prong-Grasper 4 Prong-Grasper

Catalogue No.	Grasper Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE1222-G3	3 Prong-Grasper, Sharp Hooks	2.3	230	2.8
NE1222-G4	4 Prong-Grasper, Sharp Hooks	2.3	230	2.8



Polyp Retriever

Teflon tube



3 Prong-Grasper 4 Prong-Grasper

Catalogue No.	Grasper Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NE1322-G3	3 Prong-Grasper, Blunt Eyelets	2.3	230	2.8
NE1322-G4	4 Prong-Grasper, Blunt Eyelets	2.3	230	2.8



Foreign Body Retriever

Teflon coated stainless steel spiral



Alligator Rat Tooth Crocodile

Catalogue No.	Forceps Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NEX0622-C	Alligator	2.3	180	2.8
NEX0822-C	Rat Tooth	2.3	180	2.8
NEX0822-G	Rat Tooth	2.3	230	2.8
NEX1418-A	Crocodile	1.8	120	2.0
NEX1418-C	Crocodile	1.8	180	2.0
NEX1422-C	Crocodile	2.3	180	2.8
NEX1422-G	Crocodile	2.3	230	2.8



Foreign Body Retriever Set

Teflon coated stainless steel spiral; different forceps shapes

Catalogue No.	Consists of	Forceps Shape	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
NEXSET01	3 x NEX0622-C 3 x NEX0822-C 4 x NEX1422-C	3 x Alligator 3 x Rat Tooth 4 x Crocodile	2.3	180	2.8