

SonoFlex EBUS

Endobronchial Ultrasound-Guided
Aspiration Needle – Single Use



Intuitive Locking Mechanism

For ergonomic one-hand operation

Special Surface Treatment

For superior visibility

Versions with Kink Resistant Nitinol Needle

Memory-effect prevents kinking or bending

SonoFlex EBUS TBNA – Single Use

round stylet



SonoFlex EBUS – Nitinol Needle

With atraumatic (round) Nitinol-stylet; including aspiration syringe 5–20 ml; compatible with Olympus and Fujifilm EBUS endoscopes

Catalogue No.	Description	Needle Diameter (Gauge / mm)	Outer Tube Diameter (mm)	Needle Length Adjustment (cm)
FNA-B03-2218	EBUS-guided aspiration needle including Fujifilm EBUS endoscope adapter	22/0.7	1.8	0–4

round stylet



SonoFlex EBUS – Stainless Steel Needle

With atraumatic (round) Nitinol-stylet; including aspiration syringe 5–20 ml; compatible with Olympus and Fujifilm EBUS endoscopes

Catalogue No.	Description	Needle Diameter (Gauge / mm)	Outer Tube Diameter (mm)	Needle Length Adjustment (cm)
FNA-B02-2218	EBUS-guided aspiration needle including Fujifilm EBUS endoscope adapter	22/0.7	1.8	0–4

round stylet



SonoFlex EBUS – Nitinol Needle

With atraumatic (round) Nitinol-stylet; including aspiration syringe 5–20 ml; compatible with Pentax EBUS endoscopes

Catalogue No.	Description	Needle Diameter (Gauge / mm)	Outer Tube Diameter (mm)	Needle Length Adjustment (cm)
FNA-B01-2218	EBUS-guided aspiration needle	22/0.7	1.8	0–4

LEF-GUS-SU-01-EN / LL-0390 Vers. 2.0 (11/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
 Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
 Email: sales@endo-flex.de • www.endo-flex.de