



ENDO-FLEX[®]
A MEDI-GLOBE COMPANY

Guide Wires

Single Use & Autoclavable



Teflon Coated Nitinol Wire

Kink resistance

Flexible, Hydrophilic Tip

Low friction advancement through especially challenging anatomy

Radiopaque Tip for Enhanced Visibility

Guide wire control and maneuverability under fluoroscopic imaging

Guide Wires – Single Use & Autoclavable



Nitinol Guide Wire “Hi-Flex” – Single Use

5 cm hydrophilic tip

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21725450	High Flexible Straight Tip	.025/0.64	450	5	Violet/Yellow
21735450	High Flexible Straight Tip	.035/0.89	450	5	Violet/White
21735600	High Flexible Straight Tip	.035/0.89	600	5	Violet/White



Nitinol Guide Wire “Hi-Flex” – Single Use

5 cm hydrophilic tip

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21825450	High Flexible Curved Tip	.025/0.64	450	5	Violet/Yellow
21835450	High Flexible Curved Tip	.035/0.89	450	5	Violet/White



Nitinol Guide Wire “Smart Guide” – Single Use

50 cm hydrophilic coating

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21635450	Flexible Standard Tip	.035/0.89	450	5	Violet/White



Guide Wires – Single Use & Autoclavable



Nitinol Guide Wire “Stripe Guide” – Reusable

Autoclavable; Teflon-coated

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (cm)	Colour
21525400	Standard Tip	.025/0.64	400	5	Violet/ Yellow
21535400	Standard Tip	.035/0.89	400	5	Violet/ White
21535500	Standard Tip	.035/0.89	500	5	Violet/ White
21538260	Standard Tip	.038/0.96	260	5	Violet/ White



Guide Wire Stainless Steel – Single Use

Teflon-coated

Catalogue No.	Description	Diameter (inch/mm)	Length (cm)	Flexible Part (mm)	Colour
21025400	Flexible Tip	.025/0.64	400	30	Green
21035300	Flexible Tip	.035/0.89	300	30	Green
21035400	Flexible Tip	.035/0.89	400	30	Green
21052400	Flexible Tip	.052/1.32	400	30	Green