



Multiband-Ligator

Single Use, Non Sterile



User-friendly Mounting

Adaptor-free device fixation onto endoscopes

„Shooting“ Mechanism for Rapid and Easy Operation

Facilitates band release in emergency situations

Colored Alert Band

Indicates second to last band

Multiband-Ligator – Single Use



Esophageal Multiband-Ligator MBL

Pre-loaded barrel with 6 rings (latex); additional Luer-Lock irrigation port

Catalogue No.	Suitable for Endoscope Outer Diameter (mm)	Number of Bands	Band Inner Diameter in Released Position (mm)	Length of the Introduction Catheter (cm)	Min. Accessory Channel (mm)
MBL-01-06	9.5 – 11.5	6	1.5	160	2.8

LEF-GBL-SU-00-EN / LL-0358 Vers. 1.0 (08/2020)



ENDO-FLEX GmbH

Alte Hünxer Str. 115 • 46562 Voerde • Germany
Tel.: +49 281 94400-0 • Fax: +49 281 94400-11
Email: sales@endo-flex.de • www.endo-flex.de