

# ELECTROSURGERY

## PENCILS

### DISPOSABLE PENCILS WITH PUSH-BUTTON CONTROL

ref (part n.)	cable length	Ø shaft	connector	electrode
<b>F4797</b>	320 cm	2,38mm	Valleylab	blade F4046
<b>F4797/5</b>	500 cm	2,38mm	Valleylab	blade F4046
<b>F4797/WB</b>	320 cm	2,38mm	Valleylab	none
<b>F4797/5WB</b>	500 cm	2,38mm	Valleylab	none
<b>F4797/A</b>	320 cm	2,38mm	Valleylab	needle F4048
<b>F4797/5A</b>	500 cm	2,38mm	Valleylab	needle F4048
<b>F4797/P</b>	320 cm	2,38mm	Valleylab	sphere F4044
<b>F4797/L</b>	320 cm	2,38mm	Valleylab	Blade F4050
<b>F4797/5L</b>	500 cm	2,38mm	Valleylab	Blade F4050
<b>F4797/AL</b>	320 cm	2,38mm	Valleylab	needle F4052
<b>F4797/PL</b>	320 cm	2,38mm	Valleylab	sphere F4067
<b>F4797/T</b>	320 cm	2,38mm	Valleylab	blade F4067T
<b>F4797TEF</b>	320 cm	2,38mm	Valleylab	blade F4046/TEF
<b>F4797/5TEF</b>	500 cm	2,38mm	Valleylab	blade F4046/TEF
<b>F4797/ERB</b>	320 cm	2,38mm	Erbe 5 mm	blade F4046
<b>F4797/5ERB</b>	500 cm	2,38mm	Erbe 5 mm	blade F4046
<b>F4798</b>	320 cm	2,38mm	Valleylab	blade F4046
<b>F4798/5</b>	500 cm	2,38mm	Valleylab	blade F4046



### PENCILS IN KIT

ref (part n.)	description
<b>F4798/H</b>	pencil F4797 + holster F4545
<b>F4798/5H</b>	pencils F4797/5 + holster F4545
<b>F4798/HTEF</b>	pencil F4797TEF + holster F4545
<b>F4797/75</b>	pencil F4797 + tip cleaner F7520
<b>F4797/5-75</b>	pencil F4797/5 + tip cleaner F7520
<b>F4798/75H</b>	pencil F4797 + holster F4545 + tip cleaner F7520
<b>F4798/5-75H</b>	pencil F4797/5 + tip cleaner F7520 + holster F4545
<b>F4797T/75H</b>	pencil F4797/T + holster F4545 + tip cleaner F7520
<b>F4798/5ERB-75H</b>	pencil F4797/5ERB + tip cleaner F7520 + holster F4545



F4798/75H

### DISPOSABLE PENCILS WITH SMOKE EVACUATOR SYSTEM

ref (part n.)	cable length	Ø shaft	connector	electrode
<b>F47957SE</b>	300cm	2,38 mm	Valleylab	blade F4046-F4050



F4797/SE

### DISPOSABLE PENCILS WITH ROCKER-SWITCH CONTROL

ref (part n.)	cable length	Ø shaft	connector	electrode
<b>F4795</b>	320 cm	2,38 mm	Valleylab	blade F4046
<b>F4795/5</b>	500 cm	//	//	//
<b>F4795TEF</b>	320 cm	//	//	blade F4046/TEF



F4795