

# ELECTROSURGERY

## REUSABLE HAND SWITCH PENCILS

ref (part n.)	cable length	Ø shaft	connector	n°sterilization	electrode	picture
F4141	320 cm	2,38 mm	Valleylab	30	blade F4060	
F4141/ERB	320 cm	2,38 mm	Erbe 5 mm	30	blade F4060	
F4151	500 cm	2,38 mm	Valleylab	30	blade F4060	
F4243	320 cm	2,38 mm	Valleylab	100	blade F4060	
F4244	300 cm	2,38 mm	Valleylab	100	blade F4060	
F4254	500 cm	2,38 mm	Valleylab	100	blade F4060	
F4242	300 cm	2,38 mm	Valleylab	100	blade F4060	
F4252	500 cm	2,38 mm	Valleylab	100	blade F4060	
F4246	300 cm	4 mm	Valleylab	100	blade F4060/4	
F4256	500 cm	4 mm	Valleylab	100	blade F4060/4	

## REUSABLE FOOT CONTROL PENCILS

ref (part n.)	cable length	Ø shaft	connector	n°sterilization	electrode	picture
F4814	300 cm	2,38 mm	4mm plug	100	blade F4060	

## REUSABLE ADAPTERS FOR HAND SWITCH PENCILS

ref (part n.)	pencil connector	ESU connector	compatibility	picture
F4831	Valleylab	protected 4mm plug	Berchold, Martin	
F4832	Valleylab	5mm plug	Erbe ACC/ICC	
F4833	Valleylab	4mm plug	Erbe T series	

## REUSABLE ADAPTERS FOR FOOT CONTROL PENCILS

ref (part n.)	pencil connector	ESU connector	compatibility	picture
F4830	4mm plug	8mm plug	Bovie, Soxil, Bard, Valleylab, Burdick	
F4834	4mm plug	6mm plug	Alsatom	
F4835	3,2mm plug	8mm plug	Bovie, Soxil, Bard, Danieli, Ag e Cvs, Neomed, Valleylab, Burdick, Ems	
F4836	3,2mm plug	4mm plug	Martin, Siemens, Erbe, Alsa, Laser	
F4837	3,2mm plug	4mm plug	Martin, Blendtorne, Emc, Acmicron, Erbe, Cobi	
F4838	3,2mm plug	8mm plug	Alsatom, EB SS4/MMSB	
F4839	3,2mm plug	10mm plug	Alsatom	